

MOVE ORDER REQUEST

PT. SANGHIANG PERKASA

Org CodeLAB

Move Order16865189

Move Order DescriptionMO SUBLAB MIKRO Media, Consumable, Supplement RUTIN tanggal 03122025

Transaction TypeMove Order Issue

Department

FRM-WHRPM-008.00

Need By Date03-DEC-25

UserNISA.FARHANY

Printed Date03-DEC-25

Page1

No	Item Code	Description	Qty	UOM	Sec Qty	Sec UOM	Lot	Remarks	Source Subinventry	Destination Subinventor y	COA	Timbang
1	LRRCON02	30090847;CONDU CTIVITY STANDART 1.3 US/CM (250ML);METTLER TOLEDO	2	Bot			CSKC1325J1	; RUTIN untuk kalibrasi conductivity meter	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
2	LMSLSB04	S6501-100G;STREPTOM YCIN SULFATE SALT;SIGMA	1	Bot				; RUTIN untuk supplement TOS	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
3	LMSBGA02	610009;BRILLIANT GREEN AGAR 500GR;LIOFILCHE M	4	Bot			031425502	; RUTIN untuk analisa lanjutan salmonella	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
4	LMSBPA02	CM0275B;BAIRD PARKER MEDIUM 500GR;OXOID	3	Bot			6279033	; RUTIN untuk analisa staph	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
5	LMSPDA02	CM0139B;POTAT O DEXTROSE AGAR 500 GR;OXOID	1	Bot				; RUTIN untuk analisa YM	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
6	LMSDRO02	NCM0081A;DICHL ORAN GLYCEROL (DG-18) AGAR BASE 500G,NEOGEN	1	Bot			331589	; RUTIN untuk analisa YM	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
7	LMSRVS02	610175;RAPPAPO RT VASSILIADIS BROTH 500 G;LIOFILCHEM	1	Bot			021125503	; RUTIN untuk analis lanjutan salmonella	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
8	LMSTSO02	1.05459.0500;TRY	1	Bot				; RUTIN untuk	LAB KN		03-52300-340401-	N

No	Item Code	Description	Qty	UOM	Sec Qty	Sec UOM	Lot	Remarks	Source Subinventry	Destination Subinventor y	COA	Timbang
		PTIC SOY BROTH (CASO BROTH);MERCK						analisa BC MPN			00000-1FA00-000-00000	
9	LMSLMO01	610091;DEMI FRASER BROTH 500G;LIOFILCHEM	1	Bot				; RUTIN untuk analisa listeria	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
10	LRSPOT01	POTASSIUM CHLORIDE MERCK 1.04936.1000	1	Bot			K51948436	; RUTIN untuk penyimpanan elektroda	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
11	LRRBUF03	1.09461.1000;BUF FER SOLUTION PH 9.00;1 L;MERCK	1	Bot				; RUTIN untuk kalibrasi pH meter	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
12	LANPET08	HP0006;PETRI DISH 90X15 MM;HUIDA	2	Pac				; RUTIN untuk analisa sampel liquid	LAB KN		03-52300-210501-00000-1FA00-000-00000	N
13	LANJOS01	JARUM OSE STAINLESS LOKAL	3	Pcs				; RUTIN untuk analisa lanjutan	LAB KN		03-52300-210501-00000-1FA00-000-00000	N
14	LMLEYE06	80123;EGG YOLK TELLURITE EMULSION 20% 2X50ML;LIOFILCHEM	4	Bot			010925105	; RUTIN untuk supplement BPA	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
15	LMLEYE07	80219;EGG YOLK EMULSION 2X50 ML;LIOFILCHEM	2	Bot			042325042	; RUTIN untuk supplement MYP	LAB KN		03-52300-340401-00000-1FA00-000-00000	N
16	LMSANR03	AN0035A;ANAERO GEN 3.5 L;OXOID	2	Box				; RUTIN untuk analisa bifidobacterium	LAB KN		03-52300-340401-00000-1FA00-000-00000	N

Created By

Approved By

Received By

(

)

(

)

(

)

Note To Receiver : OTH:

BK144HA; PLATE COUNT AGAR (PCA) 500GR; BIOKAR x1
BK131GC; BUFFERED PEPTON WATER 5 KG; BIOKAR x1
88029;OXIDASE TEST STICKS 50 STICK/PACK;LIOFILCHEM x2